

absorbable

MONOSWIFT™

Violet, monofilament, synthetic absorbable suture composed of Poly(glycolide-co-caprolactone) (PGCL)

comparable to Monocryl™

Benefits

- Excellent substitute for chromic gut for more predictable absorption
- Less tissue reaction than chromic gut
- Smooth tissue passage
- Excellent handling with exceptional suppleness
- Good knot security
- Non-wicking helps prevent bacterial migration along suture line
- Also available in Antibacterial PLUS line

Uses

- Short-term wound support
- Soft tissue repair/ligation
- Subcuticular closure
- GI, spay/neuter
- Not for use in cardiovascular, neurological, microsurgery or ophthalmic surgery

Wound Support

70% at 7 days
40% at 14 days
12% at 21 days

Absorption

90-120 days



The RIGHT products at the RIGHT price.



absorbable

MONOSWIFT™




Violet, monofilament, synthetic absorbable suture composed of Poly(glycolide-co-caprolactone) (PGCL)

comparable to Monocryl™

MONOSWIFT						
Needle Type	Packaging Style	Code	Size	Needle	Length	Ethicon Monocryl
PHOENIX SELECT						
▽	⬭	PS0987	0	PFSL	36"	Y987H
▽	⬭	PS0967	0	PCP-1	36"	Y267H
▽	⬭	PS0943	2/0	PFS-1	36"	Y943H
▽	⬭	PS0942	3/0	PFS-1	36"	Y942H
▽	⬭	PS0923	3/0	PFS-2	36"	Y923H
▽	⬭	PS0922	4/0	PFS-2	36"	Y922H
⊙	⬭	PS0346	0	PCT-1	36"	Y346H
⊙	⬭	PS0345	2/0	PCT-1	36"	Y345H
⊙	⬭	PS0762	2/0	PCT-2	36"	Y762H
⊙	⬭	PSL317	2/0	PSH	30"	Y317H
⊙	⬭	PS0344	3/0	PCT-1	36"	Y344H
⊙	⬭	PSL316	3/0	PSH	30"	Y316H

MONOSWIFT						
Needle Type	Packaging Style	Code	Size	Needle	Length	Ethicon Monocryl
⊙	∞	L315	4/0	NSH	30"	Y315H
⊙	∞	0310	4/0	PSH-1	30"	
▽	∞	L494	4/0	NP-3	18"	Y494G
▽	∞	L496	4/0	NPS-2	18"	Y496G
▽	∞	L493	5/0	NP-3	18"	Y493G
▽	∞	L463	5/0	NP-3	18"	Y463G
▽	∞	L495	5/0	NPS-2	18"	Y495G
ANTIBACTERIAL						
▽	∞	PL943	2/0	NFS-1	36"	MCP943H
▽	∞	PL942	3/0	NFS-1	36"	MCP942H
▽	∞	PL923	3/0	NFS-2	36"	MCP923H
▽	∞	PL922	4/0	NFS-2	36"	MCP922H
⊙	∞	PL345	2/0	NCT-1	36"	MCP333H
⊙	∞	PL316	3/0	NSH	30"	MCP316H

⊙ Taper Point ▽ Reverse Cutting ▽ Precision Reverse Cutting ∞ Suture is wound in a Figure 8 package ⬭ Suture is wound in a Racetrack package



Needle silhouettes are for display purposes only