

absorbable

Vetacryl™

Violet, coated, braided, synthetic absorbable suture composed of polyglycolic acid (PGA)

comparable to Vicryl™

Benefits

- Excellent tensile strength throughout the critical wound healing period
- Great knot security
- Exceptional handling characteristics
- Predictable absorption profile
- Coated with polycaprolactone and calcium stearate
- Also available in Antibacterial PLUS line

Uses

- General wound closure
- Soft tissue approximation and/or ligation including ophthalmic procedures
- Oral Mucosa
- Not for use in cardiovascular & neurological procedures

Wound Support

- 83% at 14 days
- 52% at 21 days
- 20% at 28 days

Absorption

60-90 days



The RIGHT products at the RIGHT price.



absorbable

Vetacryl™

Violet, coated, braided, synthetic absorbable suture composed of polyglycolic acid (PGA)

comparable to Vicryl™





VETACRYL

Needle Type	Packaging Style	Phoenix Code	Size	Needle	Length	Ethicon Vicryl
PHOENIX SELECT						
▽	⊖	PS587P	0	PFSL	30"	J587H
▽	⊖	PS453P	2/0	PFS-1	30"	J453H
▽	⊖	PS452P	3/0	PFS-1	30"	J452H
▽	⊖	PS398P	3/0	PFS-2	30"	J398H
▽	⊖	PS397P	4/0	PFS-2	30"	J397H
⊙	⊖	PS340P	0	PCT-1	30"	J340H
⊙	⊖	PS334P	0	PCT-2	30"	J334H
⊙	⊖	PS333P	2/0	PCT-2	30"	J333H
⊙	⊖	PS332P	3/0	PCT-2	30"	J332H
⊙	⊖	PS316P	3/0	PSH	30"	J316H
▽	∞	467P	0	NCP-1	30"	J467H
▽	∞	468P	1	NCP-1	30"	J468H
▽	∞	466A	2/0	NCP-1	30"	J466H
▽	∞	586P	2/0	NFSL	30"	J589H
▽	∞	443A	2/0	PFS-1	30"	J443H
▽	∞	451P	4/0	NFS-1	30"	J451H
▽	∞	391P	5/0	NFS-2	18"	J391H
▽	∞	421P	5/0	NFS-2	30"	J421H
▽	∞	492P	6/0	NP-3	18"	J492H
⊙	∞	341P	1	NCT-1	30"	J341H
⊙	∞	345P	2/0	NCT-1	36"	J345H
⊙	∞	317P	2/0	NSH	30"	J317H

VETACRYL

Needle Type	Packaging Style	Phoenix Code	Size	Needle	Length	Ethicon Vicryl
⊙	∞	339A	2/0	CT-1	30"	J339H
⊙	∞	338P	3/0	NCT-1	30"	J338H
⊙	∞	327P	3/0	NCT-3	30"	J327H
⊙	∞	310A	4/0	NSH-1	30"	J310H
ANTIBACTERIAL						
▽	∞	P453A	2/0	NFS-1	30"	VCP443H
▽	∞	P452A	3/0	NFS-1	30"	VCP452H
VISORB QUICK						
▽	∞	VQ392	4/0	NFS-2	18"	VR392
▽	∞	VQ497	3/0	PS-2	18"	VR497
▽	∞	VQ936	3/0	NPS-1	18"	VR936
▽	∞	VQ496	4/0	NPS-2	18"	VR496
▽	∞	VQ494	4/0	NP-3	18"	VR494
▽	∞	VQ493	5/0	NP-3	18"	VR493
▽	∞	VQ490	5/0	NP-1	18"	VR490
CASSETTES						
N/A	Cassette	PG10	#0	N/A	25M	
N/A	Cassette	PG01	#1	N/A	25M	
N/A	Cassette	PG02	#2	N/A	25M	
N/A	Cassette	PG03	#3	N/A	25M	
N/A	Cassette	PG20	2/0	N/A	25M	
N/A	Cassette	PG30	3/0	N/A	25M	
N/A	Cassette	PG40	4/0	N/A	25M	

⊙ Taper Point ▽ Reverse Cutting ▽ Precision Reverse Cutting ∞ Suture is wound in a Figure 8 package ⊖ Suture is wound in a Racetrack package



Needle silhouettes are for display purposes only